

Compact Click-Style Oxygen Regulators

– Designed for pediatric, homecare, hospital and emergency applications –



WESTERN

a Scott Fetzer company

Global Gas Management Technology

Enterprises | Innovator | **Medica** | Westwinds



OPA-520

Compact Click-Style Regulator with handtight CGA-540 Nut and Nipple Inlet



OPA-840

Compact Click-Style Regulator with CGA-870 Yoke Inlet

OPA-8xx | OPA-5xx

- 3000 psi Max. Inlet Pressure
- Compact aluminum body
- Click-Style flow control
- Sintered filter for additional safety and to extend regulator life
- Internal reseating relief valve protects against over-pressurization
- 1 1/2" diameter gauge
- Tapered hose barb outlet
- 2 year limited warranty

OPA SERIES COMPACT REGULATORS

Part # Flow Settings

CGA-870 Yoke

OPA-810	1/64, 1/32, 1/16, 1/8, 1/4, 1/2, 3/4, 1, 1 1/2, 2 LPM
OPA-820	1/4, 1/2, 1, 2, 3, 4, 6, 8, 10, 15 LPM
OPA-840	1/2, 1, 1 1/2, 2, 2 1/2, 3, 4, 5, 6, 8 LPM
OPA-850	1/2, 1, 2, 3, 4, 6, 8, 10, 15, 25 LPM

CGA-540 Nut and Nipple – Handtight

OPA-510	1/64, 1/32, 1/16, 1/8, 1/4, 1/2, 3/4, 1, 1 1/2, 2 LPM
OPA-520	1/4, 1/2, 1, 2, 3, 4, 6, 8, 10, 15 LPM
OPA-540	1/2, 1, 1 1/2, 2, 2 1/2, 3, 4, 5, 6, 8 LPM

⚠ WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.



Extended Click-Style Regulators



OPL-821

- 3000 psi Maximum Inlet Pressure
- Same as OPA Series PLUS
- Extended body with DISS 1240 outlet for use with a humidifier bottle

OPL SERIES EXTENDED REGULATORS

Part #	Flow Settings	Outlet Connection
CGA-870 Yoke		
OPL-821	¼, ½, 1, 2, 3, 4, 6, 8, 10, 15 LPM	DISS 1240

OPL-821

Extended Click-Style Regulator with DISS 1240 Outlet Connection Shown with a Humidifier Bottle

⚠ WARNING: This product can expose you to lead, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov.



Global Gas Management Technology

Enterprises | Innovator | **Medica** | Westwinds

Distributed by:

875 Bassett Road, Westlake, OH 44145 | US Toll Free (800) 783-7890 | (440) 871-2160 | www.westernenterprises.com

ISO 9001/ISO 13485 Registered Quality System

B20654 ECN-17925 (1806)