

PORTABLE EMERGENCY O₂ MANIFOLD

OXYGEN FOR EIGHT



©2008 Western Medica

Western's new Octi-FLO₂[™] oxygen manifold delivers portable convenience for Mass Casualty Response (MCR) Teams, EMS or triage situations. The lightweight system features our QuickClick[™] click-style flowmeters to let you quickly lock into your exact flow from 0-25 lpm for up to eight patients.

- Pressure gauge on manifold
- Standard DISS inlet and 100 psi relief valve as a safety feature
- 8 cannulas, 8 breathing masks and a compact carrying case included
- High quality craftsmanship for years of dependable service

Model MDM-8



Octi-FLO₂[™]
Portable Emergency Oxygen Manifold



**WESTERN
MEDICA**
THE LEADER IN
MEDICAL GAS EQUIPMENT

800.783.7890

INTRODUCING OCTI-FLO₂

Portable Emergency Oxygen Manifold

- Used for mass casualty, EMS, and triage situations
- Oxygen delivery for up to 8 patients
- Flowmeters are individually controlled
- 8 easy-to-read click-style flowmeters from 0-25 lpm
- Pressure gauge on manifold
- Standard DISS 1240 inlet and a 100 psi relief valve safety feature
- Lightweight and compact durable aluminum and steel construction
- Cannulas, breathing masks, and compact carrying case included
- Affordable high quality craftsmanship for years of dependable service



MDM-8

Optional Components to Complete the Emergency Oxygen Delivery System

- 0-100 psi Regulator
- DISS 20' Oxygen Hose



M1-540-PG
(0-100 psi)



C4-OXY-240-FNN-FNN
(Oxygen medical hose)

Accessories for other EMS Applications

- Wrenches
- Replacement Parts
- Washers
- T-handles
- Post Valves

MCW - 2BC



MCW - 3P



2544-2



FME101



CFM104-08
QuickClick Flowmeter



800.783.7890

Press 3 for Medical Customer Service



**WESTERN
MEDICA**
THE LEADER IN
MEDICAL GAS EQUIPMENT